POD KX 2.0

LIGAMENT KNEE BRACE

EN Instructions for use
DE Gebrauchsanweisung
FR Mode d'emploi
ES Instrucciones de uso
IT Istruzioni per l'uso
JP 使用説明書

Designed & Manufactured by /
Conçu et Fabriqué

POD Active Pty. Ltd.

Geelong, VIC 3220

AUSTRALIA

T +61 3 5244 4111
E support@podactive.com

EC REP Wellkang Ltd (www.CE-marking.eu)
29 Harley St., London W1G 9QR, UK

Made in China | Fabriqué en Chine
© Copyright 2018 POD Active | K-10506 Rev F

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ENGLISH

INSTRUCTIONS FOR USE

BEFORE USING THIS DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY.
CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE / INDICATIONS

Low to high impact/contact activities
ACL Instability
MCL/LCL Instability
Rotary and combined instabilities

Meniscus injuries
Hyperextension
Moderate ligament sprains
Prophylactic use

NOTE

The POD® knee brace is side (left and right) and size specific. For maximum comfort and support wear a knee sleeve and ensure the product fits securely.

WARNINGS AND CAUTIONS

If you experience any pain, swelling, sensation changes or any unusual reactions while using this product, consult your medical professional.

APPLICATION INFORMATION (Visit www.podactive.com for more information.)

While seated and with your knee at right angles.

1. Unfasten all four straps from the lateral/outer side of the brace.
2. Place the brace over the knee so that the hinge is located against the side of the knee joint above the mid patella line.
3. Fit **Strap #1** (directly below the knee) first as this strap is located above the calf muscle and aids in the correct positioning of the brace.
4. Fit remaining straps in numbered order.

ADJUSTMENT

Use any combination of the Hinge Pads supplied to widen or narrow the width of the brace at the knee. Control the Strap Length by adjusting the medial/Inside Velcro Strap end. The strap pads are held in position by Velcro fasteners and can be moved along the strap to increase comfort.

CARE INSTRUCTIONS

Hand wash using mild detergent and water. Air dry only.
Inspect the brace for obvious signs of wear and fatigue.

SCHEDULED SERVICE

Visit www.podactive.com for more information.

POD® HUMAN MOTION® HINGE

The patented POD® ligament is engineered with advanced materials to absorb shocks and forces on the knee.

CONSULT INSTRUCTIONS FOR USE

FOR SINGLE PATIENT USE ONLY

NOT MADE WITH
NATURAL RUBBER LATEX

⚠ WARNING ⚠

This is a supportive Knee Brace only and is not intended or guaranteed to prevent injury.
POD® Active Pty. Ltd. accepts no responsibility for any injury, or condition, which might occur while using POD® knee bracing products. If you do not agree to these conditions, contact POD® Active Pty. Ltd. before using this product.

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