

CE EC REP Wellkang Ltd
Enterprise Hub, NW Business Complex
1 Beraghmore Road, Derry, BT48 8SE
NORTHERN IRELAND

UK CA UKRP Wellkang Limited
16 Castle Street, Dover,
CT16 1PW, ENGLAND

Made in China | Fabriqué en Chine
© Copyright 2024 POD Active | K-10521 IFU REV B

EN | Instructions for use
ES | Instrucciones de uso
FR | Mode d'emploi
DE | Gebrauchsanweisung
IT | Istruzioni per l'uso
JA | 取扱説明

K8 FORGED CARBON
KNEE BRACE 3.0

pod
#1 brace in moto

CUSTOMER SERVICE:
support@podactive.com

facebook.com/podactive
twitter.com/pod_active
instagram.com/pod_active
youtube.com/podactive

podactive.com

POD Active Pty Ltd
Geelong, Victoria 3220
AUSTRALIA
T +61 3 5244 4111
E support@podactive.com

K-TYPE A
CLASS 1
MEDICAL DEVICE
EN 1621 1:2012 93/42/EEC

EN | ENGLISH

MD MEDICAL DEVICE

INTENDED USE

The device is intended for external support, stabilization and protection of the knee.

INDICATIONS FOR USE

For knee conditions that may benefit from increased anterior-posterior (AP) and medial-lateral (ML) stability and protection, such as injuries to the ACL, MCL, LCL, and PCL, and rotary and combined instabilities.
No known contraindications.

WARNINGS AND CAUTIONS

If you experience any pain, swelling, sensation changes or any unusual reactions while using this device, consult your health care professional.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

This device is not guaranteed to prevent injury.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

APPLICATION INFORMATION

While seated bend your leg to 90°.

Unfasten all four straps from the lateral/outside of the brace (**Fig. 2**).

Place the brace over the knee so that the hinge is located against the side of the knee joint above the mid patella line (**Fig. 3**).

Fit **Strap #1** (directly below the knee) first as this strap is located above the calf muscle and aids in the correct positioning of the device (**Fig. 2**).

Fit remaining straps in numbered order (**Fig. 2**).

NOTE: For maximum comfort and support wear a knee sleeve and ensure the device fits securely.

ADJUSTMENT

Control the strap length by adjusting the medial/inside hook-and-loop strap end (**Fig. 4**).

Use any combination of the **Hinge Pads** supplied to widen or narrow the width of the brace at the knee (**Fig. 5**).

The device comes standard with +0° extension stop. Install the desired **Hinge Extension Stop** by flexing the device slightly. Insert the **Hinge Extension Stop** into the hinge and push down until flush with the **Hinge Housing**, ensuring their correct orientation within the hinge. Repeat this process for the hinge on the other side, making sure to use the matching **Hinge Extension Stop** (**Fig. 6**).

The device comes with **Upper and Lower Guards** (optional) fitted, which can be easily removed from and fitted to the device with the use of an M4 hex key.

The steps for fitting and removal of the **Upper Guard** are shown (**Fig. 9**); ensure that all **Strap Loops** and fasteners are securely fitted.

The steps for fitting and removal of the **Lower Guard** are shown (**Fig. 10**); ensure that the tab and fasteners are securely fitted.

The **Elastic Cord** connects the **Impact Guards** together and ensures that they function correctly.

The steps for fitting and removal of the **Elastic Cord** are shown (**Fig. 11**); ensure that the **Hinge Pads** are correctly installed.

NOTE: When the device is used without the **Impact Guards** fitted, it is recommended that the **Impact Guard Cover** be fitted to the **Upper Frame**.

The steps for fitting and removal of the **Impact Guard Cover** are shown (**Fig. 12**); ensure that the fastener is securely fitted.

CARE INSTRUCTIONS

Hand-wash using mild detergent and rinse thoroughly.

Air dry only.

Remove foreign materials from the Hinge (e.g. dirt or grass) and clean using fresh water.

NOTE: Do not machine-wash, tumble dry, iron, bleach or wash with fabric softener.

If the device is used in salt water or chlorinated water, rinse with fresh water and air dry.

Inspect the brace for signs of wear and fatigue.

SCHEDULED SERVICE

Visit www.podactive.com for more information.

Additional information on device use, fitting, adjustment and maintenance can be found here:

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